



Universal Lateral Flow Assay Kit

For Laboratory Use

Catalog Number: AU2034-04/05

For Research Use Only. Not for use in Diagnostic Procedures.

I. Introduction

Attogene's universal lateral flow assay kit is a ready-to-use, universal test strip (dipstick), which is based on the lateral flow technology that uses gold particles containing Anti-FITC/FAM to conveniently capture nucleic acid molecules. The dipstick is designed to conveniently develop qualitative or semi-quantitative rapid test systems for detection of nucleic acids with secondary capture antibody.

Detection of nucleic acid (DNA or RNA) using this system requires the use of a biotin and fluorescein isothiocyanate or 6-carboxyfluorescein (FITC/FAM)-labelled nucleic acid for the AU2034-04 and AU2034-05-kit versions or nucleic acids containing fluorescein isothiocyanate or 6-carboxylfluorescein (FITC/FAM)-labelled digoxigenin (dig) for the AU2034-05-kit version (**Figure 1**). In creating nucleic acids using amplification, the primers should be designed so that the




Figure 1.

Concentration	Analyte Used	Control Line	Test Line	Ratio
100 pmoles	Bio-DNA-FAM	882246	1296377	1.47
1 pmoles	Bio-DNA-FAM	728786	1326008	1.82
0.981 pmoles	Bio-DNA-FAM	1248576	334428	0.27
0 pmoles	Bio-DNA-FAM	1224009	9794	0.01




Figure 2. This data demonstrates the importance of the relative amounts of molecules of a synthetic dual labeled oligo (containing biotin and FAM) with the amount of gold particles to finding optimal line intensities. GAM line captures Anti-FITC/FAM gold when Anti-FITC/FAM is not occupied with the dual labeled oligo. The more dual labeled oligo added, the binding to the Anti-FITC/FAM reaches maximal intensity at 1pmol and then diminishes due to over saturation with the dual labeled oligo above 2 pmoles, free biotinylated oligo overtakes the system to bind both the anti-FITC/FAM gold and Streptavidin capture line, while also occupying every gold particle in the system leading to a loss in intensity of both the GAM line and the Streptavidin line also called the "hook effect."

final product contains both moieties on the final product. Test line: Streptavidin, Control Line: Goat Anti-Mouse. For multiplex detection (04 kit version) of nucleic Acid (DNA or RNA) requires the use of one nucleic acid containing biotin, FITC or FAM and another containing FITC/FAM-labelled DIG: Test Line #1: Streptavidin, Line #2: anti-DIG, Line #3 Goat Anti-Mouse.

The sample containing the nucleic acids to be detected is simply mixed into the specially designed nuclease free assay running buffer in a well of the supplied 96-well plate, mixed and the dipstick is then added. Generally, the reaction is complete in 30 minutes, do not perform assay for longer than this duration. Non-specific results may occur. It is important to note that the relative stoichiometry between the nucleic acid added and the streptavidin gold is important in assay optimization (**Figure 2**). This test can be run for a ratio analysis as the more streptavidin labeled nucleic acid product that is added into the reaction the less control line is visualized. Therefore, a comparison of

the sample should be made with a strip containing only sample running buffer in each experiment. The appropriate concentration of labeled nucleic acid to use with strips is dependent upon the purity and sequence of the nucleic acid and a standard curve can be used to determine the relative ratio. A positive control dual labeled (biotin-FITC/FAM or FITC/FAM-dig) nucleic acid of your sequence that has a well characterized size, purity and concentration can be used for comparison.

Features & Benefits

- Can be used for development of a lateral flow assay for detection of a variety of different molecules such as amplified DNA products from PCR, LAMP and RPA reactions.
- No need to stripe capture antibodies
- No expensive equipment required
- Cost-effective way to screen for further downstream lateral flow assay development.

2. Kit Contents

Component Name	Volume	Storage
4.5mm Dipsticks	100 each	RT
Nucleic Acid Lateral Flow Running Buffer	20 mL	RT
Control Nucleic Acid containing biotin-FITC (AU2034-04 and AU2034-05 kit)	20 μ L	Refrigerate
Control Nucleic Acid containing FITC-DIG (AU2034-05 kit only)	20 μ L	Refrigerate
96 well plate + 4 extra wells	1 each	RT
Manual	1 each	RT

3. Storage and Stability

- The kit should be stored at 2°C - 30°C until ready to use.
- The test must remain in the sealed pouch until use.

4. Required Materials Not Supplied

- Timer - For timing use
- Centrifuge - For preparation of clear specimens
- Pipettor and pipette tips – to transfer samples and controls
- Molecule of interest containing Biotin and either FITC/FAM (04/05 version) or FITC/FAM and DIG (05 version).
- Tubes or microtiter plates to run the strips

5. Precautions

- Prior to use, ensure that the product has not expired by verifying that the date of use is prior to the expiration date on the label.
- The test strips are packaged in a foil pouch with a desiccant.
- Avoid cross-contamination of samples by using a new tube and disposable pipette tip for each sample.
- Use only Lateral Flow Kit reagents from one kit lot, as they have been adjusted in combination.
- It is good laboratory practice to use positive and negative controls to ensure proper test performance.
- Due to the hook effects, if no signal is detected in the test line, a serial dilution may be necessary to bring the nucleic acid into the appropriate concentration ratio/stoichiometry with the gold and the test line capture reagents to see the test line.

6. Procedure

Perform the following:

1. Add 150 μ L of Nucleic Acid Lateral Flow Running Buffer into a well of a 96 well plate
2. Always run a positive and negative control well with sample:
 - A. (SAMPLE) mix a designated amount (a volume 1 μ L-5 μ L are good starting points) of product into the sample running buffer. When running a LFA for the first time, we recommend trying large dilutions of sample/antigen to determine the dynamic range of the assay and keeping the volume below 15 μ L if possible.
 - B. (POSITIVE CONTROL) mix 5 μ L of the control oligonucleotide

- C. (NEGATIVE CONTROL) leave this well blank (don't add any sample or control)
3. Mix each well completely by pipetting up and down several times.
 4. Add one dipstick into each well (arrows facing up).
 5. Incubate for 30 minutes
 6. Visually analyze the strip by eye, photography or read in a lateral flow reader.

NOTE: If a higher analytical sensitivity is required, it could be helpful to increase the volume and/or concentration of the nucleic acid added into the well. Volume and concentration of analyte-specific solutions, and incubation time are always part of the individual test development.

NOTE: Control oligonucleotide line should yield a test and control line signal on the strip within 15 minutes.

7. Interpretation of Results

This test is a lateral flow assay containing test lines that are dependent on the concentration of biotin-FITC/FAM (for AU2034-04 and AU2034-05 kit) and/or FITC/FAM-DIG (for AU2034-05 kit) labeled nucleic acids in the sample.

What to expect at the test lines:

The higher the concentration of nucleic acid in the sample the higher the intensity of the test line compared to the strip lacking nucleic acid (negative control strip).

What to expect at the control line:

The intensity of the control will decrease as the test line increases.

Determinations made using strips which have dried for more or less than the required time may be inaccurate, as line intensities may vary with drying time.

8. Additional Analysis

If necessary, positive samples can be confirmed by analyzing using a nucleic acid analysis technique such as agarose or acrylamide gels. A lateral flow reader may also be employed to generate numerical readings.