



Blank 4.5mm Lateral Flow Strips for LFA Development

Catalog Number: **AU2054**

Product Data Sheet

Introduction:

Attogene's Blank 4.5mm Lateral Flow Strips for LFA Development is a convenient, simple, and ready to be used out of the box product designed to help labs with limited facilities engage in proof-of-concept LFA development. We provide a product that facilitates equal, if not superior results to leading brands products at a significantly improved cost. With a high barrier to entry, we recognize that it is important to provide laboratories with accessible components that will allow for more innovation in this sector.

Kit Contents:

1. 100 4.5X60mm strips
2. Two strip cannisters with an integrated desiccant.
3. One 96 well plate

Component	Description	Properties
CO-UTCP80-17	Untreated Glass Fiber Conjugate Pad	Thickness: 0.42mm Wicking Rate: 2.2 (s/2 cm)
CO-WM110-22	Cellulose Fiber Wick	Thickness: 1.10mm
CO-SS40-60	Nitrocellulose Membrane	Mylar Thickness: 100µm Membrane Thickness: 100µm 4cm Wicking Time: 90 sec. Protein Binding: Highest protein binding capacity available

Usage guide:

- 1. Capture antibody/antigen conjugate application.**
 - a. Spot reagents between 0.1-2mg/mL. There may be instances where going beyond this range may yield improved sensitivity, however, this is a typical concentration range.
 - b. Use a protein such as BSA, or a saccharide such as sucrose to stabilize reagents on the nitrocellulose.
 - c. Keep spotting volumes around 0.5-1uL total.
 - d. Depending on the application, it is recommended to dry the strips at 2-8°C for a period of 12-24 hours or at 37°C for anywhere between 1-12 hours.



2. Running the direct LFA:

- a. For usage directly in a microwell, sample/running buffer amounts between 80-200uL is advised.
- b. Mix your sample and running buffer with 1-10µL of 10 OD of the appropriate gold conjugate. For usage with a fluorescent reporter particle such as Quantum Dots, Europium Beads, or Phycoerythrin proper amounts must be determined via titration.
- c. Insert the strip's conjugate pad into the microwell.
- d. Typical LFA running times range from 5-30 minutes. This must be determined through experimentation.

Related Attogene Products Catalog Number:

1. Protein labeling and purification

- a. **CO2027** Microspin Desalting Columns 7K MWCO
- b. **CO2028** Protein Biotinylation Kit
- c. **AU2059** Lateral Flow Starter Kit

2. Consumables

- a. **CO2022** Lateral Flow Cassettes
- b. **CO2033** 96-well Plate
- c. **AU2035** RNase Free Buffer Screening Pack
- d. **AU2054** 4.5mm X 60mm Strips for Lateral Flow Development

3. Reporter Reagents

- a. **AU2012** 10nm Colloidal Gold
- b. **AU2013** 30nm Colloidal Gold
- c. **AU2014** 40nm Colloidal Gold
- d. **AU2015** 60nm Colloidal Gold
- e. **AU2016** Protein A Colloidal Gold for Lateral Flow
- f. **AU2017** Streptavidin Colloidal Gold for Lateral Flow
- g. **AU2043** Streptavidin CdSe/ZnS Quantum Dots for Lateral Flow
- h. **AU2050** Streptavidin Europium Beads for Lateral Flow
- i. **AU2055** Streptavidin Phycoerythrin for Lateral Flow

Troubleshooting:

If it becomes evident that spotted reagents aren't adhering to the nitrocellulose surface. It may be beneficial to add methanol or ethanol at a very low concentration to improve binding kinetics. Ensure reagent's biocompatibility before adding any organic solvents. Protein or molecules may also need to be attached to a carrier protein such as BSA.

Non-specific results may be induced by a variety of factors such as the LFA running time was too long, buffer is incompatible, reagents are not compatible, or reagents were not blocked sufficiently/correctly etc.

Storage:

All components should be stored at temperatures ranging from 2°C-RT. To ensure longevity of all components, avoid exposure to extreme temperatures and leave the strips in their provided canisters which contain an integrated desiccant component.

Product Safety and Handling:

Components provided in this kit are for R&D purposes only. This kit is not for use in diagnostic procedures.