



Rapid Test Components Plastic Cassettes


Product Lateral Flow Cassettes – Customization – Perfect for LFA Product
Intended Use Development With Standard Card Stock
Product Variant Development/Manufacture of Lateral Flow Immunoassay Test Kit
 CO2022-02

Property	All Type
Appearance	White
Fitment	OK

Type	Description
CO2022-02	5 mm plastic cassette

Specifications

REPRESENTATIVE DIAGRAM OF PLASTIC CASSETTE - NOT TO SCALE



Parameters of Plastic Cassettes	CO2022-02 Dimensions
A. Sample Hole Diameter	4.1 mm
B. Width of Strip	5.0 mm
C. Length of Strip	60 mm
D1. Distance from sample hole to Pressure point (sample pad side)	10 mm
D2. Distance from sample hole to Pressure point (Absorbent pad side)	35.5 mm
D3. Distance between sample hole and lower side of square hole	10.7 mm
L1. Width of the device	20.5
L2. Length of device	70.6
Length of window	16.8
Width of window	4.2
NC Membrane length on the strip	25

Storage Condition Store and handle in accordance with all current regulations and standards. Store in original packing at room temperature. Avoid direct sunlight.
Caution: Material should be kept away from heat radiators.